



Syntron Bioresearch, Inc.

QuikScreen® Multi X Drug Cup Test

ONE STEP ONSITE DRUG CUP

Catalog # 60XXX

CLLAWAIVED

— Instructions —



INTENDED USE

The **QuikScreen® Multi X Drug Cup test** is an immunochromatographic assay for rapid, qualitative detection of drug combinations and their principal metabolites in urine at specified cut-off concentrations. In the **QuikScreen® Multi X Drug Cup test**, X may denote any number of drugs. These drug combinations may be composed from any of the following drugs, at the noted cut-off concentrations:

DRUG CLASS	ABBREVIATIONS	SENSITIVITY
AMPHETAMINE	AMP	1000 ng/ml
BARBITURATES	BAR	300 ng/ml
BENZODIAZEPINES	BZD	300 ng/ml
COCAINE	COC	300 ng/ml
MARIJUANA	THC	50 ng/ml
METHADONE	MAD	300 ng/ml
METHAMPHETAMINE	MET	1000 ng/ml
OPIATES	OPI	2000 ng/ml
OXYCODONE	OXY	100 ng/ml
PHENCYCLIDINE	PCP	25 ng/ml
TRICYCLIC ANTIDEPRESSANT (2)	TCA	1000 ng/ml

Note: The test provides only preliminary data which should be confirmed by other methods such as gas chromatography/mass spectrometry (GC/MS). Clinical considerations and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

SUMMARY AND EXPLANATION OF THE TEST

The **QuikScreen® Multi X Drug Cup test** is an easy, fast, qualitative, visually read competitive binding immunoassay method for screening without the need of instrumentation. The method employs unique mixture of antibodies to selectively identify the drugs of abuse and their metabolites in test samples with a high degree of sensitivity.

Drug abuse remains a growing social and economical concern in many developed and developing countries throughout the world. The above stated drugs are among the most frequently abused illicit drugs, according to the U.S. Substance Abuse and Mental Health Services Administration. Opiates are among a class of heavily abused prescription drugs.

The sensitivity of the **QuikScreen® Multi X Drug Cup test** is set as required for the screening immunoassays of these drugs in the reference guidelines set by the U.S. Substance Abuse and Mental Health Services Administration (SAMHSA) and the U.S. Department of Health and Human Services.

PRINCIPLE OF THE TEST

The **QuikScreen® Multi X Drug Cup test** is a competitive binding immunoassay in which drug and drug metabolites in a urine sample compete with immobilized drug conjugate for limited labeled antibody binding sites. By utilizing antibodies that are specific to different drug classes, the test permits independent, simultaneous detection of any of the drug combinations from a single sample. The approximate run time is 5 minutes.

In the assay procedure, urine mixes with labeled antibody-dye conjugate and migrates along a porous membrane. When the concentration of a given drug is below the detection limit of the test, unbound antibody-dye conjugate binds to antigen conjugate immobilized on the membrane, producing a rose-pink color band in the appropriate Test Zone for that drug. Conversely, when the drug level is at or above the detection limit, free drug competes with the immobilized antigen conjugate on the membrane by binding to antibody-dye conjugate, forming an antigen-antibody complex, preventing the development of a rose-pink color band.

Regardless of the drug levels in the sample, a rose pink-color band is produced in each Control Zone (top bands) by a parallel immunochemical reaction. These bands serve as built-in quality



control measures by demonstrating antibody recognition, verifying that the reagents are chemically active.

REAGENTS AND MATERIAL PROVIDED

1. Test Devices	Contains dye-conjugated antibody and immobilized antigen in protein matrix with sodium azide.
2. Test Instructions	Catalog# PI-60XXX-CLIA
Optional:	
3. Negative Control I	Contains buffered protein solution with sodium azide. Cat.# 4010N
4. Amphetamine Positive Control	Cat.# 11120P
5. Barbiturates Positive Control	Cat.# 18040P
6. Benzodiazepines Positive Control	Cat.# 18020P
7. Cocaine Positive Control	Cat.# 12000P
8. Marijuana Positive Control	Cat.# 13020P
9. Methadone Positive Control	Cat.# 19020P
10. Methamphetamine Positive Control	Cat.# 11320P
11. Opiates Positive Control	Cat.# 11220P
12. Oxycodone Positive Control	Cat.# 19080P
13. Phencyclidine Positive Control	Cat.# 14020P
14. Tricyclic Antidepressant Positive Control	Cat.# 19092P

MATERIALS REQUIRED BUT NOT PROVIDED

1. Clock or timer.
2. Specimen collection containers.

WARNINGS AND PRECAUTIONS

1. For *in-vitro* diagnostic use.
2. Do not use the test device beyond the expiration date.
3. Collect urine specimen directly into the test cup. Ensure that the sample amount meets the minimum level as indicated on the side of the test cup.
4. Read the results at 5 minutes. Do not interpret results after 30 minutes.

STORAGE AND STABILITY

Store test kit below 28°C; **do not freeze**. If stored at 2°-8°C, allow the test kit to reach room temperature (15°-28°C) before performing the test. Refer to the expiration date for stability.

SPECIMEN COLLECTION AND PREPARATION

Fresh urine specimens should be collected directly into the cup. The **QuikScreen® Multi X Drug Cup test** device employs a **thermal strip which should be checked immediately** after collection to validate urine specimen. SAMHSA regulations specify that any temperature below 90.5° F must be considered adulterated. No additives or preservatives are required.

Note: Urine specimens can be transferred from a urine collection container into **QuikScreen® Multi X Drug Cup test cup**, if necessary.

TEST PROCEDURE

1. Do not break the seal of the pouch until ready to begin testing.
2. Remove the test cup from the foil pouch.
3. Collect urine specimen directly into the test cup. Ensure that the sample amount meets the minimum level as indicated on the side of the test cup.
4. Read the results at 5 minutes. Do not interpret results after 30 minutes.

Note: The result must be interpreted at five minutes. Waiting more than five minutes may cause the reading to be inaccurate. To avoid confusion, discard the test device after interpreting the result.

INTERPRETATION OF RESULTS

Positive: A *rose-pink* band is visible in each control zone (top band). No color band appearing in the appropriate test zone (bottom band) indicates a preliminary positive result for the corresponding drug of that specific test zone. Send urine specimen to a certified laboratory for confirmation.

Negative: A *rose-pink* band is visible in each control zone and the appropriate test zone, indicating that the concentration of the corresponding drug of that specific test zone is below the detection limit of the test.

Invalid: If a color band is not visible in each of the control zones, the test is invalid. Another test should be run to re-evaluate the specimen.


For Professional Use



Syntron Bioresearch, Inc. • 2774 Loker Avenue West • Carlsbad, CA 92010 • USA
Emergo Europe • Molenstraat 15 2513 BH • The Hague • Netherlands +31.70.345.8570

Doc. ID: PI-60XXX-CLIA
Revision: 02a
Eff. Date: 09/30/11
DC03773

Note: There is no meaning attributed to line color intensity or width.



QUALITY CONTROL

An internal procedure control has been incorporated into the test to insure proper kit performance and reliability.

The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS OF THE TEST

1. This product is designed to be used for the detection of drugs of abuse and their metabolites in human urine only.
2. Although the test is very accurate, there is the possibility false results will occur due to the presence of interfering substances in the specimen sample.
3. The test is a qualitative screening assay and is not suggested for quantitative determination of drug levels in urine, or the level of intoxication.
4. Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, can cause erroneous test results regardless of the analysis method used.
5. If adulteration is suspected, obtain another urine specimen.

PERFORMANCE CHARACTERISTICS

1. **Sensitivity.** The QuikScreen® Multi X Drug Cup test detects drugs of abuse and their major metabolites in urine at concentrations equal to or greater than the cut-off level for the specific drug, which is suggested by the U.S. Substance Abuse and Mental Health Services Administration (SAMHSA) for the immunoassay method.
2. **Specificity.** A study was conducted with the QuikScreen® Multi X Drug Cup test to determine the cross-reactivity of drug-related compounds with the test. Substances listed in **Table I** produced results approximately equivalent to the cutoff levels. A separate study was conducted to determine the cross-reactivity of non-related compounds with the test at concentrations much higher than normally found in the urine of people using or abusing them. No cross reactivity was detected with the substances listed in **Table II**.

Table I: Concentrations of drug-related compounds showing positive response approximately equivalent to the cut-off set for the test:

The following Amphetamine-related substances yield positive results for Amphetamine at 1000 ng/ml cut-off level:

d-Amphetamine	1000 ng/ml
l-Amphetamine	25,000 ng/ml
Δ-l-Amphetamine	10,000 ng/ml
β-Phenylethylamine	180,000 ng/ml
Thyramine	100,000 ng/ml
(±) 3,4-Methylenedioxyamphetamine-HCl (MDA)	1200 ng/ml

The following Barbiturate-related substances yield a positive result for Barbiturates at 300 ng/ml cut-off level:

Allobarbitol	600 ng/ml
Amobarbitol	600 ng/ml
Barbitol	300 ng/ml
Butobarbitol	300 ng/ml
Butalbital	300 ng/ml
Pentobarbitol	300 ng/ml
Phenobarbitol	300 ng/ml
Secobarbitol	300 ng/ml

The following Benzodiazepine-related substances yield positive results for Benzodiazepines at 300 ng/ml cut-off level:

Alprazolam	600 ng/ml
Bromazepam	100 ng/ml
Chlordiazepoxide	300 ng/ml
Clobazam	300 ng/ml
Clonazepam	300 ng/ml
Clorazepate	200 ng/ml
Delorazepam	3,000 ng/ml

Diazepam	300 ng/ml
Estazolam	300 ng/ml
Flunitrazepam	300 ng/ml
Flurazepam	150 ng/ml
Lorazepam	500 ng/ml
Lormetazepam	500 ng/ml
Nitrazepam	250 ng/ml
Nordiazepam	150 ng/ml
Oxazepam	300 ng/ml
Prazepam	1,500 ng/ml
Temazepam	150 ng/ml
Triazolam	200 ng/ml

The following Cocaine-related substances yield positive results for Cocaine at 300 ng/ml cut-off level:

Benzoylcegonine	300 ng/ml
Cocaine	300 ng/ml

The following Marijuana-related substances yield positive results for Marijuana at 50 ng/ml cut-off level:

Cannabinol	10,000 ng/ml
11-nor-Δ-8-THC-9-COOH	50 ng/ml
11-nor-Δ-9-THC-9-COOH	50 ng/ml
Δ-8-THC	7500 ng/ml
Δ-9-THC	10,000 ng/ml
11-hydroxy-Δ-9-THC	2500 ng/ml

The following Methadone-related substances yield positive results for Methadone at 300 ng/ml cut-off level:

Methadone	300 ng/ml
Doxylamine	50,000 ng/ml
EDDP (2 Ethylidene-1,5-dimethyl 1-3,3-Diphenylpyrrolidin)	100,000 ng/ml
Methadol	25,000 ng/ml
Perphenazine	75,000 ng/ml
Protriptyline	2,000 ng/ml
Trimipramine	10,000 ng/ml

The following Methamphetamine-related substances yield positive results for Methamphetamine at 1000 ng/ml cut-off level:

(+) Methamphetamine	1000 ng/ml
(±)3,4Methylenedioxyamphetamine (MDMA)	1000 ng/ml
(±)3,4Methylenedioxyamphetamine (MDA)	200,000 ng/ml
d-Amphetamine Sulfate	200,000 ng/ml
l-Amphetamine Sulfate	200,000 ng/ml
Δ-l-Amphetamine Sulfate	200,000 ng/ml

The following Opiates-related substances yield a positive result for Opiates at 2000 ng/ml cut-off level:

Morphine	2000 ng/ml
Morphine Sulfate Pentahydrate	2000 ng/ml
Morphine-3-β-D Glucuronide	2000 ng/ml
Codeine	2000 ng/ml
Heroin	2000 ng/ml
Levorphanol	4000 ng/ml
Ranitidine	100,000 ng/ml
6-Acetylmorphine	50 ng/ml

The following Oxycodone-related substances yield positive results for Oxycodone at 100 ng/ml cut-off level:

Oxycodone-HCl	100 ng/ml
Codeine	700 ng/ml
Hydrocodone	500 ng/ml
Hydromorphone	1,500 ng/ml
Morphine-Sulfate	7,000 ng/ml
Morphine-3-b-D-Glucuronide	40,000 ng/ml
Norcodeine	40,000 ng/ml
Oxymorphone	300 ng/ml

The following Phencyclidine-related substances yield a positive result for Phencyclidine at 25 ng/ml cut-off level:

Phencyclidine	25 ng/ml
Tenocyclidine	2000 ng/ml

The following Tricyclic Antidepressant-related substances yield positive results for Tricyclic Antidepressant at 1000 ng/ml cut-off level:

Amitriptyline	1,000 ng/ml
Cyclobenzaprine	1,500 ng/ml
Clomipramine	5,000 ng/ml
Desipramine	600 ng/ml
Doxepin	1,000 ng/ml
Imipramine	600 ng/ml
Notriptyline	1,000 ng/ml
Nordoxepin	1,000 ng/ml

Table II: Compounds tested and found not to cross-react with the test at a 100 µg / ml concentrate in urine.

Acetaminophen	Ecgonine methyl ester	Oxazepam
Acetone	EDDP	Omega-3-Fatty Acid
Acetyl Salicylic Acid	+ Ephedrine	Penicillin G
albumin	- Ephedrine	Pentobarbital
Amikacin	(+/-) Epinephrine	Penphenazine
Amitriptyline	Erythromycin	Perphenazine
Amobarbital	Ethanol	Phenalzine
Amphetamine	Fentanyl	Phencyclidine
D Amphetamine	Fluxetine	pheniramine
Ampicillin	Furosemide	Phenobarbital
Arterenal	Gentisic acid	Phentermine
l-Ascorbic Acid (Vitamin C)	Glucosamine	l-Phenylephrine
Aspartame	Glucose	Phenylethylamine-a
Aspirin	Guaiacol Glyceryl Ether	(+/-)-Phenylpropanolamine
Atropine	DL-Hematropine	Primidone
Benzocaine	Hemoglobin	Procaine
Benzoic Acid	Histamine	Promethazine
Benzoylcegonine	Homatrophine	d-Propoxyphene
Benzoylcegonine HCl	Hydrochlorothiazide	Pseudoephedrine
(+)- Brompheniramine	Hydrocodone	Quinine
Buprenorphine	Hydromophone	Quinine antidine
buprenorphine-3-β-d-glucuronide	Ibuprofen	Quinidine
Butabartital	Imipramine	Salicylic Acid
Caffeine	Isoproterenol	Secobarbital
Camphor	Ketamine	Sodium Chloride
(+)-Chlorpheniramine	Lidocaine	Sulindac
(+/-)-Chlorpheniramine	Maprotiline	Sustiva
Chlorprothixene	3,4±MDA	Tenocyclidine
Chloroquine	(+/-)MDMA	Tetracycline
Chlorpromazine	Meperidine	Delta-9- Tetrahydrocannabinol
Cocaine	Methadone	Tetrahydrozoline
Codeine	Methamphetamine	Theophylline
Cortisone	Methanol	Thioridazine
(-)-Cotinine	Methaqualone	11-nor-Δ-THC-9-COOH (10 ug/ml)
Creatinine	Methylphenidate	11-nor-Δ-THC-9-COOH (10 ug/ml)
r-Cyclodextrin	Morphine	Tramadol
Cyclobenzaprine	Morphine-3-β-D-Glucuronide	d(+)-Trehalose
Deoxyephedrine	(1S,2S)-(-)-N-methyl-ephedrine	Trifluoperazine
Dextromethorphan	Naloxone	Trimethobenzamide
Diazepam	Naltrexone	Tyramine
Digitoxin	b-Naphthaleneacetic acid	Triprolidine Hydrochloride
Digoxin	(+/-) Naproxen	
4-Dimethylaminoantipyrine	Neomycin	
Diphenhydramine	Niacinamide	
5,5-Diphenylhydantoin	Nicotene	
Dopamine	Nicotinic Acid	
Doxylamine	Nor-Buprenorphine	
Ecgonine	Noscapine Hydrochloride	
Ecgonine HCl	Oxalic Acid	

In order to examine potential naturally occurring interfering substances normally contained in urine. Drug free urine and drug positive urine were spiked with various potential interfering substances. Both samples were tested with QuikScreen® Multi X Drug Screen Test Cup Device. No cross-reaction was noted by any of the following substances at the concentrations list in the following table.

Table III - Natural Occurring Compounds in Urine and the Effect on QuikScreen® Multi X Drug Cup test

Analyte	Range	Effect	
		Positive*	Negative**
Ascorbic	300 mg/dl		
Bilirubin	1.0 mg/dl	None	None
Creatine	500 mg/dl	None	None
Glucose	1500 mg/dl	None	None
Hemoglobin	300 mg/dl	None	None
Potassium	110 mEq/dl	None	None
Human Serum Albumin	500 mg/dl	None	None
Globulin	1500 mg/dl	None	None
Sodium chloride	6000 mg/dl	None	None
Uric Acid	23 mg/dl	None	None
Cholesterol	500 mg/dl	None	None

*Concentration of Positive Drug Control = Amphetamine 1250 ng/ml, Methamphetamine 1250 ng/ml, Opiates 2500 ng/ml, Cocaine 375ng/ml, THC 63 ng/ml, Phencyclidine (PCP) 32 ng/ml, Benzodiazepine (450 ng/ml), Barbiturate (450 ng/ml), Methadone (450 ng/ml), TCA (1250 ng/ml), Oxycodone (100 ng/ml), Propoxyphene (300 ng/ml).

** Concentration of Drug [Drug Free urine] = 0 ng/ml

Effects of prolonged specimen exposure to the test device: In order to determine if there were any significant affects on the specimen by prolonged exposure to the test device, a study on the **QuikScreen® Multi X Drug Cup test** was performed using in-house urine control with GC/MS value assignment. The test specimens were subjected to a time zero (0) GC/MS evaluation. The test specimens were then applied to the **QuikScreen® Multi X Drug Cup test** such that the fluid level was midway between urine level marks and moderately shaken for a period of 10 minutes. The **QuikScreen® Multi X Drug Cup test** with the test specimens were stored for 50 hours at room temperature(15-30°) C.

Samples for GC/MS analysis were taken at times 0, 12, 36 and 60 hours. Statistically there was no significant change in the concentrations reported for any of the analytes at any time period. Based upon the GC/MS data, it may safe to conclude that there were no significant changes in the analyte concentrations of specimens that could be related to the device or the test strips contained in the device.

Accuracy: The accuracy of the **QuikScreen® Multi X Drug Cup test** was tested in a clinical trial of urine samples submitted to a SAMHSA certified laboratory. The laboratory used Syva® EMIT II as their screening procedure. All positive samples by either screening method were confirmed by GC/MS. The relative sensitivity results by either GCMS is summarized as follows:

3.1 AMPHETAMINE (AMP) 1000 NG/ML CUT-OFF LEVEL		
	GC/MS Positive	GC/MS Negative
QuikScreen® Positive	47	2
QuikScreen® Negative	0	61

When compared to GC/Mass the relative sensitivity was computed to be 47/47 or 100%. The relative specificity was computed to be 61/63 or 97%. The concordance of the combined data with respect to GC/Mass was 108/110 or 98%.

3.2 BARBITURATES (BAR) 300 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	39	2
QuikScreen® Negative	0	52

When compared to GC/Mass the relative sensitivity was computed to be 39/39 or 100%. The relative specificity was computed to be 52/54 or 96%. The concordance of the combined data with respect to GC/Mass was 91/93 or 98%.

3.3 BENZODIAZEPINE (BZD) 300NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	42	2
QuikScreen® Negative	0	52

When compared to GC/Mass the relative sensitivity was computed to be 42/42 or 100%. The relative specificity was computed to be 52/54 or 96%. The concordance of the combined data with respect to GC/Mass was 94/96 or 98%.

3.4 COCAINE (COC) 300 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	42	1
QuikScreen® Negative	0	53

When compared to GC/Mass the relative sensitivity was computed to be 42/42 or 100%. The relative specificity was computed to be 53/54 or 98%. The concordance of the combined data with respect to GC/Mass was 95/96 or 99%.

3.5 MARIJUANA (THC) 50 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	42	1
QuikScreen® Negative	0	53

When compared to GC/Mass the relative sensitivity was computed to be 42/42 or 100%. The relative specificity was computed to be 53/54 or 98%. The concordance of the combined data with respect to GC/Mass was 95/96 or 99%.

3.6 METHADONE (MAD) 300 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	42	0
QuikScreen® Negative	0	52

When compared to GC/Mass the relative sensitivity was computed to be 42/42 or 100%. The relative specificity was computed to be 52/52 or 100%. The concordance of the combined data with respect to GC/Mass was 94/94 or 100%.

3.7 METHAMPHETAMINE (MET) 1000 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	41	1
QuikScreen® Negative	0	59

When compared to GC/Mass the relative sensitivity was computed to be 41/41 or 100%. The relative specificity was computed to be 59/60 or 98%. The concordance of the combined data with respect to GC/Mass was 100/101 or 99%.

3.8 OPIATES (OPI) 2000 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	41	1
QuikScreen® Negative	0	54

When compared to GC/Mass the relative sensitivity was computed to be 41/41 or 100%. The relative specificity was computed to be 54/55 or 98.2%. The concordance of the combined data with respect to GC/Mass was 95/96 or 99%.

3.9 OXYCODONE (OXY) 100 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	55	1
QuikScreen® Negative	0	43

When compared to GC/Mass the relative sensitivity was computed to be 55/55 or 100%. The relative specificity was computed to be 43/44 or 98%. The concordance of the combined data with respect to GC/Mass was 98/99 or 99%.

3.10 PHENCYCLIDINE (PCP) 25 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	40	4
QuikScreen® Negative	0	53

When compared to GC/Mass the relative sensitivity was computed to be 40/40 or 100%. The relative specificity was computed to be 53/57 or 93%. The concordance of the combined data with respect to GC/Mass was 93/97 or 96%.

3.11 TRICYCLIC ANTIDEPRESSANT (TCA) 1000 NG/ML CUT-OFF LEVEL

	<u>GC/MS Positive</u>	<u>GC/MS Negative</u>
QuikScreen® Positive	41	1
QuikScreen® Negative	0	42

When compared to GC/Mass the relative sensitivity was computed to be 41/41 or 100%. The relative specificity was computed to be 42/43 or 98%. The concordance of the combined data with respect to GC/Mass was 83/84 or 99%.

BIBLIOGRAPHY**GENERAL**

- Baselt, R.C., Disposition of Toxic Drugs and Chemicals in Man, 2nd Ed., *Biomedical Publ.*, Davis, CA, p.488 (1982).
- Cody, J.T., Schwarzhoff, R., *J. Anal. Toxicol.*, 17: 2630 (1993).
- Urine Testing for Drugs of Abuse, *NIDA Research Monograph* 73, (1986).
- Dasguspta, A., Saldana, S., Kinnaman, G., Smith M., Johansen, K., *Clin. Chem.*, 39(1):104-108 (1993).
- Department of Health and Human Services, *Fed. Regist.*, 53(69): 11970-11989 (1988), (1989).
- FDA Guidance for Labeling Urine Drugs of Abuse Screening Testing*, Kshitij Mohan, 7/21/1987.
- Blum, K., *Handbook of Abusable Drugs*, Gardner Press, Inc., New York, 1st Ed., (1984).
- Tietz, N.W.: *Clinical Guide to Laboratory Tests*; W.B. Saunders Company, (1976).

AMPHETAMINE/METHAMPHETAMINE TEST

- Ellerbe, P., Long, T., Welch, M.J., *J. Anal. Toxicol.*, 17: 165-170 (1993).

BARBITURATE TEST

- Hoffmann, F.E. *A Handbook of Drug and Alcohol Abuse: The Biomedical Aspects.* Oxford University Press, New York, 1983.
- Wyngarrden, J.B., Smith, L.H. (eds) *Cecil Textbook of Medicine.* W.B. Saunders Co., Philadelphia, 1988, pp. 53-54.
- Ellenhorn, M.J., Barceloux, D.J. *Medical Toxicology.* Elsevier Science Publishing Co., New York, 1988, pp. 575-580.
- Gorodetzky, Z.V. *Detection of Drugs of Abuse in Biological Fluids.* In: Martin, W.R. (ed.) *Drug Addiction.* I. Springer-Verlag, New York, 1977, pp. 319-409.

BENZODIAZEPINES TEST

- S.J., Mule and G.A. Gasella. "Quantitation and Confirmation of the Diazolo and Triazolobenzodiazepines in Human Urine by GC/MS," *J. Anal. Toxicol.*, Vol. 13 (1989), pp. 179-184.
- C. Drouet-Coassolo, C. Aubert, P. Coassolo, and J. Cano. "Capillary GC/MS Method for the Identification and qualification of some Benzodiazepines and Their Unconjugated Metabolites in Plasma," *J. Chromatogr.*, Vol. 487 (1989), pp. 295-311.

- C. Ballanto, V. Reggi, G. Tognoni, et al., "Benzodiazepines: Clinical Pharmacology and Therapeutic Use," *Drugs*, Vol. 19 (1980): 195-219.
- A.J. Giannini and E.A. Slaby. *Drugs of Abuse* (MEDEC Books: 1989), Ch. 4, pp. 59-80.
- M.A. Pear and L. Kopjak, "The Screening and Quantitation of Diazepam, Flurazepam, Chlordiazepoxide, and Their Metabolites in Blood and Plasma by Electron-Capture Gas Chromatography and High Pressure Liquid Chromatography," *J. Forensic Sci.*, Vol. 24(1979) pp. 46-54.
- H. Schultz, "Modern Screening Strategies in Analytical Toxicology with Special Regard to New Benzodiazepines," *J. Legal. Med.*, Vol. 100 (1988) pp. 19-37.
- J.G. Langner, B.K. Gan, R.H. Liu, L.D. Baugh, P. Chand, J.L. Weng, et al. "Enzymatic Digestion, Solid Phase Extraction, and GC/MS of Derivatized Intact Oxazepam in Urine," *Clin. Chem.*, Vol. 37 (1991) pp.1596-1601.
- R.L. Fitzgerald, D.A. Rexin, and D.A. Herold. "Benzodiazepine Analysis by Negative Chemical Ionization GC/MS," *J. Anal. Toxicol.*, Vol. 17, pp.342-347.

COCAINE TEST

- D.W. Hoyt et al., *J. Am. Med. Assoc.*, Vol. 258 (1987), pp. 504-509.
- R.R. MacGregor, J.S. Fowler, and A.P. Wolf. *J. of Chromatography*, Vol. 590 (1992), pp. 354-58.
- E.J. Cone, D. Yousefnejad, and S.L. Dickerson. *J. of Forensic Sciences*, Vol. 35, No.4 (1990), pp. 786-91.
- E.J. Cone et al., *J. of Forensic Sciences*, Vol. 34, No. 1 (1989), pp. 15-31.
- B. Holmstedt and A. Fredga. *J. of Ethnopharmacology*, Vol. 3 (1981), pp. 113-47.
- T. Inaba and J. Can. *Physiol. Pharmacol.*, Vol. 67 (1989), pp. 1154-57
- A.R. Jeffcoat et al., *Drug Metabolism and Disposition*, Vol. 17, No. 2 (1989), pp. 153-59.

MARIJUANA TEST

- Johansson, E., Gillespie, H.K., Halldin, M.M., *J. Anal. Toxicol.*, 14: 176-180 (1990).
- El Sohly, M.A., Jones, A.B., El Sohly, H.N., *J. Anal. Toxicol.*, 14: 277-279 (1990); 14: 227-279 (1990).
- Foltz, R.L., Sunshine, I., *J. Anal. Toxicol.*, 14: 375-378 (1990).
- Wimbish, G.H., Johnson, K.D., *J. Anal. Toxicol.*, 14: 292-295 (1990).
- Nakamura, G.R., Meeks, R.D., Stall, W.J., *J. Forensic Sci.*, 35(4): 792-796 (1990).
- Jenkins, A.J., Mills, L.C., Darwin, W.D., Huestis, M.A., Cone, E.J., Mitchell, J.M., *J. Anal. Toxicol.*, 17:6.
- Hollister, L.E., Kanter, S.L., Board, R.D., Green D.E. *Res. Com. Chem. Pathol. Pharmacol*, 8: 579-584 (1974).

METHADONE TEST

- B.J. Rounsaville, M.M. Weissman, P.H. Rosenberger, et al. "Diagnosis and Symptoms of Depression in Opiate Addicts: Course and Relationship to Treatment Outcome," *Arch. Gen. Psychiatry* 39 (1982): 151-56.
- B.J. Rounsaville and H.D. Kleber, "Psychiatric Disorders and the Course of Diagnostic Stability," in S.M. Mirin, e.d., *Substance Abuse and Psychopathology* (American Psychiatric Press, Washington, D.C., 1984): 134-51.
- C.A. Dackis and M.S. Gold. "Depression in Opiate Addicts," *ibid*, pp. 20-40.
- J.B. Wyngarrden, LH Jr., eds., *Cecil Textbook of Medicine* (WB Saunders Co.: Philadelphia, 1988), pp. 57 and 60.
- M.J. Ellenhorn and D.G. Barceloux. *Medical Toxicology* (Elsevier Science Publishing Co., Inc., NY, 1988), pp. 714-18.

OPIATE TEST

- Huang, W., Andollo, W., Hearn W.L., *J. Anal. Toxicol.*, 16: 307-310 (1992).
- Cone, E.J., Dickerson, S., Paul, B.D., Mitchell, J.M., *J. Anal. Toxicol.*, 17: 156-164 (1993).
- Glare, P.A., Walsh, T.D., and Pippenger, C.E., *Ther. Drug Monit.*, 13: 226-232 (1991).
- Walsh, T.D., Cheater, F.M., *Pharm. J.*, 10: 525-527 (1983).
- Mitchell, J.M., Paul, B.D., Welch, P., Cone, E.J., *J. Anal. Toxicol.*, 15: 49-53 (1991).

OXYCODONE TEST

- Bolan, E.A., et al., Synergy between opioid ligands: evidence for functional interactions among opioid receptor subtypes. *J. Pharmacol. Exp. Ther.* 303, 557-562 (2002).
- Ross, F.B., and Smith, M.T., The intrinsic antinociceptive effects of oxycodone appear to opioid receptor mediated. *Pain* 73, 151-157 (1997).

PHENCYCLIDINE TEST

- Froelich, P.E., Gross, G.: Separation and Detection of Phencyclidine in Urine by Gas Chromatography, *J. Chromatograph* 1977; 137; 135-143.
- Gupta R.C. et al: Determination of Phencyclidine in Urine and Illicit Street Drug Samples, *Clin. Toxicol* 1975; 8:611-621.
- Oellerich, M.: Enzyme Immunoassays In Clinical Chemistry: Present Status and Trends, *J. Clin. Chem. Biochem.*, 1980;18:197-208.
- Anilineo, Pitts, F.N.: Phencyclidine (PCP): A Review and Perspectives. *CRC Crit. Rev. Toxicol* 1982; 10:145-177.

TRICYCLIC ANTIDEPRESSANT TEST

- Drug Facts and Comparisons, 55th edition, *St. Louis, A. Wolters Kluwer Co.* 2001., 902-911.
- Wong R., The Effect of adulterants on urine screen for drug of abuse: Detection by an on-site dipstick Device, *Am. Clin. Lab.*, 2002: 21(3); 14-18.
- Fed. Register. Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs: 53, 69, 11970-11979 (1988).
- U.S. Department of Transportation, Drug Testing Procedures Handbook.

