

ORAL-VIEW™ Saliva Multi-Drug of Abuse Test

5 DRUG PROCEDURE CARD

1

Remove test device from pouch. Write ID on back of device then uncap to expose collection pad.

result window

cap

2

result window

collection pad

Insert collection pad into mouth with test device slanted downwards. Hold edge of plastic firmly between the lips with ample distance from result window. Do not bite, chew or pull on collection pad.

3

Leave in mouth for 1-3 minutes or until pink color begins to appear in result window. Remove device from mouth immediately.

Insert no more than 1 cm of plastic between subject's lips

4

control line

REPLACE CAP.
Place device on flat surface.
Start timing when C (control) line appears

5-7 min

5

Read results at 5-7 minutes after the C line first appears. Do NOT read after 7 minutes

<p>C</p> <p>THC</p> <p>OPI</p> <p>COC</p> <p>BZD</p> <p>AMP</p> <p>Negative Colored line(s) appear(s) next to the respective drugs</p>	<p>C</p> <p>THC</p> <p>OPI</p> <p>COC ← Positive for COC</p> <p>BZD</p> <p>AMP</p> <p>Positive No colored lines appear next to the respective drugs</p>	<p>C</p> <p>THC</p> <p>OPI</p> <p>COC</p> <p>BZD</p> <p>AMP</p> <p>Invalid No C line appears</p>
---	--	---

ORAL-VIEW™ Saliva Multi-Drug of Abuse Test

5 DRUG PROCEDURE CARD

1

Remove test device from pouch. Write ID on back of device then uncap to expose collection pad.

result window

cap

2

result window

collection pad

Insert collection pad into mouth with test device slanted downwards. Hold edge of plastic firmly between the lips with ample distance from result window. Do not bite, chew or pull on collection pad.

3

Leave in mouth for 1-3 minutes or until pink color begins to appear in result window. Remove device from mouth immediately.

Insert no more than 1 cm of plastic between subject's lips

4

control line

REPLACE CAP.
Place device on flat surface.
Start timing when C (control) line appears

5-7 min

5

Read results at 5-7 minutes after the C line first appears. Do NOT read after 7 minutes

<p>C</p> <p>THC</p> <p>OPI</p> <p>COC</p> <p>BZD</p> <p>AMP</p> <p>Negative Colored line(s) appear(s) next to the respective drugs</p>	<p>C</p> <p>THC</p> <p>OPI</p> <p>COC ← Positive for COC</p> <p>BZD</p> <p>AMP</p> <p>Positive No colored lines appear next to the respective drugs</p>	<p>C</p> <p>THC</p> <p>OPI</p> <p>COC</p> <p>BZD</p> <p>AMP</p> <p>Invalid No C line appears</p>
---	--	---